

We put science to work.™

# TechBriefs

## Savannah River National Laboratory

U.S. DEPARTMENT OF ENERGY • SAVANNAH RIVER SITE • AIKEN • SC

[srnl.doe.gov](http://srnl.doe.gov)

### Benefits

- > More sensitive than conventional or high-resolution-UT-only systems
- > Shortens test duration, reducing large-scale test costs (simulant, operations, etc.)
- > Minimizes process flow stream disruption
- > Provides a correlation between UT and mass loss data
- > Provides quantifiable data on erosion patterns, scar depth and pit morphology

### Applications and Industries

- > Corrosion coupon providers (Metals Samples)
- > Corrosion Test Laboratories (DNV GL Labs and National Labs)
- > End users such as oil and gas industry and chemical industry

### Intellectual Property

- > US Patent SRNS-99 (SRS-16-013 and 12-027)
- > Testing underway at Florida International University (supporting Hanford Tank Farm)
- > Custom as-required

### Contact Information

E-mail: [partnerships@srnl.doe.gov](mailto:partnerships@srnl.doe.gov)



SRNL-L1100-2018-00009

## Replaceable Ultrasonic Thickness/Mass Loss Measurement Coupon


Savannah River National Laboratory has developed a replaceable coupon that provides both high-resolution ultrasonic thickness (UT) and mass loss measurements of erosion and corrosion in critical locations throughout a process system without disrupting the process flow stream.

### Description

Incorporating mass loss as well as high-resolution UT measurement capability in a single device provides more sensitive measures of wear than conventional systems or systems providing high-resolution UT measurements alone. In addition, the unique size and design allow installation in vulnerable locations in a process system to provide timely information about the condition of the actively operating system. Because it can be used in an actively operating system, this device can dramatically reduce large-scale testing costs.

The coupon can provide in situ high-resolution UT measurements in most locations (piping, tanks, ducts, injection points, etc.) in a process system, including those with different pipe sizes/wall thicknesses. The design is optimized to allow accurate placement of a high-resolution UT pencil probe, which provides repeatable and more precise real-time wear data. This innovative design provides a means to evaluate critical evidence of corrosion and erosion without significantly disrupting the process flow stream.

The coupon can also be removed periodically to nondestructively evaluate mass loss as well as surface features resulting from corrosion (general and localized) and erosion. Once the coupon is removed, mass loss can be verified by independent techniques, such as ultrasonic and conventional micrometer measurement. An optical microscope or scanning electron microscope can be used to investigate active corrosion mechanisms, such as pitting and crevice corrosion. Finally, scar depth and orientation can be accurately evaluated using a surface profilometer or a laser confocal microscope.



Schematic showing coupon installed in a pipe



Removable Ultrasonic Thickness/Mass Loss Coupon



UT probe inserted into a coupon installed in the extrados of a 3 inch, schedule 40, 90 degree stainless elbow